

ULK3-unc-51-like kinase 3 (C. elegans) Gene

Reference: PTXBC056423

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	serine/threonine-protein kinase ULK3; unc-51 like kinase 3		
For research use only.			

GENE DESCRIPTION	unc-51-like kinase 3 (C. elegans)
GENE ID	25989
ORIGIN SPECIES	Human
ACCESSIONS	BC056423
VECTOR	pDONR223
NCBI SYMBOL	ULK3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcggggccggctgggtccccgcgctggacgcttcacacacgagcgcctgggagcggcagctacgccacgggtgtacaaggcctacgccaagaaggacactcgtgaagtggtagccataaagtgttagccaagaaaagtctgaacaaggcatcgggtgagaacct
cctcacggagattgagatcctcaaggcattcgacatccccacattgtgcagctgaagactttcagtgaggacgtgacaatatctacctcatcattggagttttgcgcaggggcgacgtgtctcgtctcatccatacccgaggaltctgcctgagaagggtggcggtgtctcatgcagcaa
ttagctagcgcctgcaattcctgcatgaacggaatatctctacctggatctgaagccacagaacattctactgagctccttggagaagccccacctaataactggcagactttgtttcgcacaacacatgtccccgtggatgagaagcacgtgtctccgtggctccccctctacatggccc
ccgagatgggtgcccagcggcagatgacgcccgtggacctctgggtccatgggggtcatctgtatgtgtgagacacttttccctgtttcaccctga