

## ZMAT2-zinc finger, matrin type 2 Gene

Reference: PTXBC056668

Product type: DNA Storage Condition: -20°C  
Size/Quantity: 2 ug Delivery Condition: +4°C, Blue Ice  
Synonyms: Ptg-12; Snu23; zinc finger matrin-type protein 2; zinc finger matrin-type 2  
**For research use only.**

---

<b>GENE DESCRIPTION</b>	zinc finger, matrin type 2
<b>GENE ID</b>	153527
<b>ORIGIN SPECIES</b>	Human
<b>ACCESSIONS</b>	BC056668
<b>VECTOR</b>	pDONR223
<b>NCBI SYMBOL</b>	ZMAT2
<b>LENGTH</b>	ProteoGenix bp
<b>SEQUENCE PRIMERS</b>	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

### ORF (Only for reference, product provided without STOP codon)

```
atggcgtgggcagcgggacaaaaactggacttgcggaaaagtggaacaagatgaatgagaaactgcgcgagaagaggctcacggaagagagagaaaaagatgaaaaccagtcagcctgtcaagcagagctttacggcataggactacaaggaggacttga
atccaagcttgggaagacaattgtcattaccaagacaaccctcaatctgagatgggaggatattactgcaatgtctgtgactgtgtggaaggactccatcaactttctggatcacattaatgaaaagaacatcagagaaacctggcatgtctatgcgttggaaactccacctggatc
aggtgaagaacgttttgaggtcaacaagaagaagtggaagagaagcagaaggattatgatttgaggaaaggatgaaggagctcagagaagaggagaaaaggccaaagcgtacaaagaagagaacagaaggagaagaaaaggaggctgagggactgacattgagg
aggacgatgagatggcagctgtgatggccttctctgcttggcttccaccaagaagattactga
```