

DCUN1D2-DCN1, defective in cullin neddylation 1, domain containing 2 (S. cerevisiae)
Gene

Reference: PTXBC056669

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C13orf17; DCN1-like protein 2; DCN1, defective in cullin neddylation 1, domain containing 2; DCUN1 domain-containing protein 2; defective in cullin neddylation protein 1-like protein 2; defective in cullin neddylation 1 domain containing 2		

For research use only.

GENE DESCRIPTION	DCN1, defective in cullin neddylation 1, domain containing 2 (S. cerevisiae)
GENE ID	55208
ORIGIN SPECIES	Human
ACCESSIONS	BC056669
VECTOR	pDONR223
NCBI SYMBOL	DCUN1D2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgcataagcttaaatcgtctcagaaggacaaggtccgcccagttatggcgtgcaactcaggtggcgagagaactgctactactgcttaacgcagaatgagtgagactagacgagggcacggacagcttctccaaaccagactcgctccacagggagtcctatgcggaacgtgtg
gacaagaagaagctggagcggctgtacggcaggtacaaagatccacaagatgaaacaaaattggagtcgagtgattcaacagttttgatgatctgagcctggatcctccagatcagtgattggtcatagcgtggaagttcagggcagcaactcagtgtaatttagcagaagg
aatttctagatggcatgacagaacttgggtgtgacagcatggagaagctaaaggtcttctgccaagactggagcagggagctgaaggacacagccaagttaagattttatcagtttaccttcaccttcgtaagaaccagggcagaagggttagacttagaaatggctgttgcgtattgg
aaattagtttatctggaaggtttaatttttagatctctggaacacattcttaattggaatcacaaaagatcaattccaaggacacctggaacctctgctggacttggaaacatgattgcggatgatgtcttaactacgatgaagaaggagcttggccgttcttatagatgattttagaat
atgcacggccagtagtcacaggtggaacgcagcctttctag