

CAPS2-calciphosine 2 Gene

Reference: PTXBC056672

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms:	UG0636c06; calciphosin-2; calciphosphine 2; calciphosine 2		

For research use only.

GENE DESCRIPTION	calciphosine 2
GENE ID	84698
ORIGIN SPECIES	Human
ACCESSIONS	BC056672
VECTOR	pDONR223
NCBI SYMBOL	CAPS2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

```
atgataaccacttatctagggtctcatcagtgatccagagcaaaattagccattgagcaaaaagaagtgatcatatcctccagattcaaatgacacctctcatttagaaaaagaactacatgaaacaaagataagaactcattacatcaactgaaatgtcctttcataaatt
acagttgatggtagatgatcacgaacaaatgctccctttatcaaaaagcatttatagtcacgtgtggacgaagaaaagaaaacaataccgacttggatgatttatgttggtgcaacctgacatggtgagttctgatcatctcagcctccagaagcatcaaaagaaacacatta
ctaaactccgaatcacaatatgatcaaatagcttggatctctcaaaactcctctatggaacaggagatgatataatcattcaagaaccaatgataggctggtctcctcaagcaattcaagatgtgctaaaagaaaaactacataaaagaggtgtctgattttgactgattgggaaat
ttcaacagttggacaaggaaatggacttttagataagcagatftaagcaagctctaaaagtgttccattagaaggtctgaaaaggatfttagtctgcatggctaatctgaatgacaatggcaatggcaaggtgattatggagaattcaaacgtggtatttgggaaatgaatga
atacagaaatcatatgtcgaaggcctttatgaaactggattcaacaaaagtggcagtgctgctattataacataagaaaatgttactgtgcaagaagcattctcaagtaattcaggttga
```