

C16orf87-chromosome 16 open reading frame 87 Gene

Reference: PTXBC056676

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	UPF0547 protein C16orf87; chromosome 16 open reading frame 87		

For research use only.

GENE DESCRIPTION	chromosome 16 open reading frame 87
GENE ID	388272
ORIGIN SPECIES	Human
ACCESSIONS	BC056676
VECTOR	pDONR223
NCBI SYMBOL	C16orf87
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtctgcaactcgagccaagaagtgaaagtgccaccaaatcatgccccgagtcgacccaacaggttcctgtgcatgtaaatcatgtcctgtggttacatatttattagcagaaaactttaaatgcaaaactcagagaaatcacccttctacagaaaacaagcatgaggccaaga
ggaggcgaacagagagagtaggagagagaagataaattctacagtaataaagatttagaaaacagaaagaggtctcgaagtaacagccattcagatcatatcagacgaggaagaggaagacctaaaagtcacatctgccccaaaaacatgaggaagaaagagagaacagaaaaag
gaaattgacatctatgctaactgtctgatgaaaaggcttctgttttcagtcgccttggcagaaataatagaaaaattatcaatcaaaacttattctctga