

DNAJC2-DnaJ (Hsp40) homolog, subfamily C, member 2 Gene

Reference: PTXBC056682

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	MPHOSPH11; MPP11; ZRF1; ZUO1; dnaJ homolog subfamily C member 2; DnaJ (Hsp40) homolog, subfamily C, member 2; M-phase phosphoprotein 11; zuotin-related factor 1; DnaJ heat shock protein family (Hsp40) member C2		

For research use only.

GENE DESCRIPTION	DnaJ (Hsp40) homolog, subfamily C, member 2
GENE ID	27000
ORIGIN SPECIES	Human
ACCESSIONS	BC056682
VECTOR	pDONR223
NCBI SYMBOL	DNAJC2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgctgcttctgccaagcgcgcggacggccggggcaccgccatcacccacgctctgacctctgacctctacactctgtcaagttgaacctgtgggaagatggtttgaagctttgttaagaggagaacagaaatgcttctgctcttttcaggaaactggaggataaagaggtatccgag
gaatcagaagatgaagaattgcagttggaagatttccatgctgaaaacacttgatcccaagactggagaacacagatcattatgcagttcttgacttggccatgtgagatacaaggctacacagagacagatcaaaagcagctcataaagcaafggttttaaacatcacccagacaaa
cggaaagcagctggtgaaccaataaagaaggagataatgactacttctgataactaaaggatccccaggtgattcttatcagaagctcatcctacagtttaattcatcttatgaatgttatctgatccagtga