

HBQ1-hemoglobin, theta 1 Gene

Reference: PTXBC056686

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	HBQ; hemoglobin subunit theta-1; hemoglobin theta-1 chain; hemoglobin, theta 1; theta-1-globin; hemoglobin subunit theta 1		

For research use only.

GENE DESCRIPTION	hemoglobin, theta 1
GENE ID	3049
ORIGIN SPECIES	Human
ACCESSIONS	BC056686
VECTOR	pDONR223
NCBI SYMBOL	HBQ1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcgctgtccgggagaccggcgctgtgctgcgcctgtggaagaagctggcgagcaacgtcggcgtctacacgacagggccctggaaaggaccttctggtttccccgccacgaagacctacttctccacctggacctgagccccggctcctcacaagtcagagccac
ggccagaaggtggcgacgctgagcctcgcctggagcgcctggacgacctacccacgcgctgtccgcgctgagccacctgcacgcgtgccagctgcgagtgagccggccagcttccagctcctggccactgctgctgtaacctcgcgggactacccggagactt
cagccccgcgtgcaggcgtcgtggacaagttcctgagccacgttatctcggcgtgtgttccgagtagcctga