

YWHAQ-tyrosine 3-monooxygenase/tryptophan 5-monooxygenase activation protein, theta polypeptide Gene

Reference: PTXBC056867

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	1C5; HS1; 14-3-3 protein theta; 14-3-3 protein T-cell; 14-3-3 protein tau; protein, theta; tyrosine 3-monooxygenase/tryptophan 5-monooxygenase activation protein, theta polypeptide; tyrosine 3-monooxygenase/tryptophan 5-monooxygenase activation protein th		

For research use only.

GENE DESCRIPTION	tyrosine 3-monooxygenase/tryptophan 5-monooxygenase activation protein, theta polypeptide
GENE ID	10971
ORIGIN SPECIES	Human
ACCESSIONS	BC056867
VECTOR	pDONR223
NCBI SYMBOL	YWHAQ
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggagaagactgagctgatccagaaggccaagctggccgagcagccgagcgtctacgacacatggccacctgcatgaaggcagtgaccgagcagggcccgagctgtccacgaggagcgcaacctgctctcgtggcctacaagaacgtggtcggggcccgaggtccgctggagggctcatctctagcatcgagcagaagaccgacacctccgacaagaagttgcagctgattaaggactatcgggagaaagtgagtcggagctgagatccatctgcaccacggtgctggaattgttggataaatatftaagccaatgcaactaatccagagagtaaggctctctatctgaaaatgaagggtgattacttccggtaccttgcgaagttgctgtggtgatcgaacaaacgataataatcccaaggagcttaccagagggcatttgatataagcaagaagagatgcaacccacacaccaatccgctggggcttgccttaactttctgtattttactatgagattcttaataaccagagcttgcctgcacgctgcttaaacgctttgatgagccattgctgaacttgatacactgaatgaagactatacaaaagacagaccctcatcatgcagttgcttagagacaacctaactttggacatcagacagtgacaggagaagaatgtgatgcccagagaaggctgaaaactaa