

C7orf44-chromosome 7 open reading frame 44 Gene

Reference: PTXBC056884

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C7orf44; MITRAC15; cytochrome c oxidase assembly factor 1 homolog; cytochrome c oxidase assembly protein 1 homolog; mitochondrial translation regulation assembly intermediate of cytochrome c oxidase protein of 15 kDa		

For research use only.

GENE DESCRIPTION	chromosome 7 open reading frame 44
GENE ID	55744
ORIGIN SPECIES	Human
ACCESSIONS	BC056884
VECTOR	pDONR223
NCBI SYMBOL	C7orf44
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgatgtggcaaaagtatgcaggaagcaggcggcgaatgcctctgggagcaaggatcctttccacgggtgtgttctatgccggggccttgcattgtgtattacctcattcaaaagttcattcaggcgtttatattacaagttggcagtgaggcagctgcagagccatcccaggcacagga
agcctctggccctctctcaacatcattatcaagctcatcgacagggaacacgttgacattgtgatgccaagtgaagattcctgtctctggatccaaatcagaggccttctctacgtccactatccagagggtggcccttcagagggtggcaccttgacgaggtcttttagagctc
aaggatggtcagcagattcctgtgtcaagctcagtgagggaacgggtgatgaagtgaaggaaggagtag