

LOC51233-hypothetical protein LOC51233 Gene

Reference: PTXBC056888

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C22orf43; KB-208E9.1; aspartate-rich protein 1; aspartate rich 1		
For research use only.			

GENE DESCRIPTION	hypothetical protein LOC51233
GENE ID	51233
ORIGIN SPECIES	Human
ACCESSIONS	BC056888
VECTOR	pDONR223
NCBI SYMBOL	LOC51233
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggggaatatactgacctgtgtatcaactcccactgtggctggccagggggaaggacgaccctgttatgaatctgatactgatatattatgagactgtggctgtgcaacatcagaatccactactgtagagcctggcaagctggatgtgggagccacggaggccaagacctgcagca
catcagcaacccaaagatgcccacaggtcccccctgaggaccgctgagtttaaaatttctgcatcaagtggaggaagacaatgatgatgccaagattttaccatcacctgtccagggtttcttgaggacaacctgagtttagtatgctaccacgaagtgaagatgatgactgtgatgat
gatgatgatgatgccagattttaccgtcacgtgccagggtggctgttaccgggtttgatagcagttctgttctcttgaggacaacctgagtttagtgctaccacgaagtgaagatgatgactgtgatgatgatgatgatgccagattttaccgtcacgtgtccaggctgtgtctgaagat
agcctgttttaagatgctcactgagacacaaagatgaagaagaagaagatgatgatccacataacagctcggatagaagtgactgacgtggagagcttaagtgatgaagagattcatccagggtga