

ARPC5-actin related protein 2/3 complex, subunit 5, 16kDa Gene

Reference: PTXBC057237

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	ARC16; dJ127C7.3; p16-Arc; actin-related protein 2/3 complex subunit 5; Arp2/3 protein complex subunit p16; actin related protein 2/3 complex, subunit 5, 16kDa; arp2/3 complex 16 kDa subunit; actin related protein 2/3 complex subunit 5		

For research use only.

GENE DESCRIPTION	actin related protein 2/3 complex, subunit 5, 16kDa
GENE ID	10092
ORIGIN SPECIES	Human
ACCESSIONS	BC057237
VECTOR	pDONR223
NCBI SYMBOL	ARPC5
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtcgaagaacacagtgctcggcccgcttcggaagggtggacgtggatgaatgacgagaacaagttcgtggacgaagaagatggggcgacggccaggccggcccgacgaggcgaggtggattcctgctcggccattccatcacaggaaacatgacagctgcctaca
ggcagctctgaagaacccctatcaacaccaagagtcaggcagtgaggaccggcgagcagcattgtcttgaaggctcctcttttaaagctaataatagataagaaaggcagttcaatctctggacaagaatggtgtggatctcctaatgaagtatatataaaggattgagagcccg
tctgacaatagcagtgctatgtactgcaatggcatgaaaaggcacttgcctgctggaggagtagggctcattgtctgtcttgactgcaagaaaactgtgtag