

NINJ2-ninjurin 2 Gene

Reference: PTXBC057766

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	ninjurin-2; nerve injury-induced protein 2; ninjurin 2		
For research use only.			

GENE DESCRIPTION	ninjurin 2
GENE ID	4815
ORIGIN SPECIES	Human
ACCESSIONS	BC057766
VECTOR	pDONR223
NCBI SYMBOL	NINJ2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggaatcagcaagagaaacatcgacctcaacctggaagctcgcagcccgagcccatcaacctgaaccattacgccaccaagaagcgtggcggagagcatgctggacgtggccctgttcatgtccaacgccatgcggctgaaggcgggtgctggagcagggaccatc
ctctcaactactacaccacctgtgtcaccctcatcagcctctctctgtctctgcaggtgggtcatcgggtgctctgctcgtggtcattgcacggctgaacctgaatgaggtagaaaagcagtggtgactcaaccagctcaacaacgcagccaccatctgtctcttctcactgtgtcatcaatgttttc
attacagcttcggggcacataaaacaggggttctggtgcaaggccctcaagggaatcctctctga