

FAM26D-family with sequence similarity 26, member D Gene

Reference: PTXBC057769

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	protein FAM26D; C6orf78; family with sequence similarity 26 member D		
For research use only.			

GENE DESCRIPTION	family with sequence similarity 26, member D
GENE ID	221301
ORIGIN SPECIES	Human
ACCESSIONS	BC057769
VECTOR	pDONR223
NCBI SYMBOL	FAM26D
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgctgggttgatttgcacacattggcaaccattgctgcttagtctctgctgtgtggcaagtgtgctctccctcactctctgcaacattgctactggaccagccacctccagaatgagagagaactctttgaacaagcagcagagcagcactctcggtcctcatgatgcacgcata
aagaagctatttggcttcatcccgaggagtgagacgtcaaacacatccgcatctctctgtcaggactggaagatattcagtaccactctttatgcatgggtgatgacttgaaggtcactatagcttcttggaaatagggtggatgaggataatgaggagacagatcaagaggtatt
gaattaaaaccttga