



Antibodies and Proteins

Espace Européen de l'Entreprise

67300 Schiltigheim, France

3 90 20 54 70 - Fax +33 (0)3

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

GPR78-G protein-coupled receptor 78 Gene

Reference: PTXBC057778

DNA

Storage Condition:

-20°C

2 ug

Delivery Condition:

+4°C, Blue Ice

G-protein coupled receptor 78; G protein-coupled receptor 78

For research use only.

GENE DESCRIPTION	G protein-coupled receptor 78
GENE ID	27201
ORIGIN SPECIES	Human
ACCESSIONS	BC057778
VECTOR	pDONR223
NCBI SYMBOL	GPR78
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggccccccggcaggcgctgctgctggcggttctctctgtgatgtgtaacgctggcgtgctatcaacgactgtgtctgtttgttgccgtacagcgtgactgcgcactcgagcctcaggcgtctctctgtgtaaatctgtctctggcccaactgctgctggcgccgctggacatgc
cttccacgtgctgctgctgtatgcctggcggtggcgccggcgccggcgctcaagaaatcattgttctctggacacctcttcgctggcgtcaacgctggcgctgagcgtggcggtgctggcggtgctggcggtgctggcggtcttccactgctgacgtgcgtacgcgcggcagctgcgcac
cgcgatgctcggcgctgctgctggcggtgtctcgtagggacagctgctgtcttccatcagcgtgcagctctctgttgcgtgacagacgcttctcgtctgttgcgtgctgctgcgcgcggcagctgagcgtgctgcgtgcgcgttcacgcttccacgctcagcgtctgcgcac
gtctgctgctgcgcgtcgctgggtgctgctcctcatcctgcacagctggcgccgggtgctgcacgacactgcgcagcgcattgacacacgctcatgaagcgtcgtcgtgctgctgcgcacgtcaccctcagcgtgctgcgcacgtgcgcacgagcgtgcgcacgcgc
gcaccagagaagatggcatctgctatfctgcacacttctcatctgttccccgtatgcatgacagcagctggcggtgctgcttgcctctgacctgcacgcggcagcgtggcgctcatctcagaagctgctgcactcatcagaaaggcgtgctgcacgcgttcacgtgcgcg
gctgttcgcgaagtcttggcgcatggtgcacggcgctgctgaagagaaccccgcgccgacatccacccatgacagactgcttggatgtggcgccgcatggtgcaccagctgctgaagagaaccccgccggcagctgcacccacacagcgtctgttgacacagagaatgattcgtg
ctctcagcagacacactga

More info: Please contact us for any question: contact@proteogenix.fr

Page 1/1