



Reference: PTXBC057780

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	RKST1; SGLT6; SMIT2; sodium/myo-inositol cotransporter 2; Na(+)/myo-inositol cotransporter 2; homolog of rabbit KST1; sodium-dependent glucose cotransporter; sodium/glucose cotransporter KST1; sodium/myo-inositol transporter 2; solute carrier family 5 (so		

For research use only.

GENE DESCRIPTION	solute carrier family 5 (sodium/glucose cotransporter), member 11
GENE ID	115584
ORIGIN SPECIES	Human
ACCESSIONS	BC057780
VECTOR	pDONR223
NCBI SYMBOL	SLC5A11
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]