



Reference: PTXBC057781

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	SERKAL; WNT-4; protein Wnt-4; wingless-type MMTV integration site family, member 4; Wnt family member 4		

For research use only.

GENE DESCRIPTION	wingless-type MMTV integration site family, member 4
GENE ID	54361
ORIGIN SPECIES	Human
ACCESSIONS	BC057781
VECTOR	pDONR223
NCBI SYMBOL	WNT4
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]

Page 1/1