



Reference: PTXBC057786

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	SLRR5B; podocan-like protein 1; podocan like 1		
For research use only.			

GENE DESCRIPTION	podocan-like 1
GENE ID	79883
ORIGIN SPECIES	Human
ACCESSIONS	BC057786
VECTOR	pDONR223
NCBI SYMBOL	PODNL1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]