

TMEM180-transmembrane protein 180 Gene

Reference: PTXBC057787

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	TMEM180; C10orf77; bA18I14.8; transmembrane protein 180; major facilitator superfamily domain containing 13A		

For research use only.

GENE DESCRIPTION	transmembrane protein 180
GENE ID	79847
ORIGIN SPECIES	Human
ACCESSIONS	BC057787
VECTOR	pDONR223
NCBI SYMBOL	TMEM180
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggacctggtgcaggcttccactgccactcaacagcaactcttccctcttctcggagcatctgtgtccgaccatattcccttccacgggtccatctgttgggctctctatgtcgtcccatctcaaacctctacttctgtccctgtgccggcgtggggcgtctacgcggtg
gtgcgggggctcttctctcaagctgggacttagcctgctcatgttgtgcccggcccgaccacctcagcctgctgtcccttcatgcccagcaaccgcttctcactgaggccacctgtaagctgctgaccttgggtgctactgacctggttagacgagacctggtgctgaaccaccgc
aagcaggcagcctcggcactcctcttggcatggttgccttgggtgaccaagccaggccagaccttggcccgctgctgggacactggctgctctgtttctacacaggtcatgacctcttccagcagctccctcataacccctgtggggagtgcccatccctggccagagccccagctccagc
ccctgcacaggccccgacgctccgagggtgcttctacctgctggtgctgtgcccacacctgtgctctgtgcagctcttccctgggtccagtcacgctgcatgggagacgctgcacatggtcaaggccagcgcagaaacctgtcacaggcccaaacctggatgttaagat
ggtgtga