

LYPD5-LY6/PLAUR domain containing 5 Gene

Reference: PTXBC057816

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	PRO4356; ly6/PLAUR domain-containing protein 5; metastasis-associated protein; LY6/PLAUR domain containing 5		

For research use only.

GENE DESCRIPTION	LY6/PLAUR domain containing 5
GENE ID	284348
ORIGIN SPECIES	Human
ACCESSIONS	BC057816
VECTOR	pDONR223
NCBI SYMBOL	LYPD5
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgggggtccccagagtcattctgctctcttggggctgcactctgacctgacagggccaagccctgcagtgctacagctttgacacacactactttggccctttgacctcaggccatgaagctgccagcatctctgtctcatgagtgctttgaggctatctgtcttgacacc
gggtatcgcgccggtgacctggtgaggaaggctgctggaccgggctctgctgggcccagacgcaatgaacgagcgctgcccagactactcggtggtgctgctgcaacaactgacaaatgaacgccacctcatgactcatgacgcccctcccaacctgagcc
aagcaccgacccgccgacgctcagcggcgccgagtgctacgctgtatcggggtccaccaggatgactgcgctatcggcaggtcccagcaggtccagtgctcaccaggaccagaccgctgctccaggccaatggcagaatgacagttggcaattctcagtcctgtgtacatcag
aacctgccaccggccctctgaccaccgagggcagcaccagccctggacagccatcgaacctcagggtctctgctgtgaggggtacctctgcaacaggaaatccatgaccagcccttcaccagtgcttcagccaccacctccccgagcactacaggtcctgcccctgctctc
ccagtcctctgctggtgggctctcagcatag