

ANKRD13A-ankyrin repeat domain 13A Gene

Reference: PTXBC057817

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	ANKRD13; NY-REN-25; ankyrin repeat domain-containing protein 13A; NY-REN-25 antigen; ankyrin repeat domain 13A		

For research use only.

GENE DESCRIPTION	ankyrin repeat domain 13A
GENE ID	88455
ORIGIN SPECIES	Human
ACCESSIONS	BC057817
VECTOR	pDONR223
NCBI SYMBOL	ANKRD13A
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtctcgcgcctgcgacgcccgcgaccactaccctgcacctcctagtctggaaaaacgactaccggcagctcgagaaggagctgcagggccgaatgtggagctgtggaccacgaggtcgaaattattgcattctgtttctctgggacatttgaatctgctcgagttact
ccgacataaagcagatgtgacaaaagaaatccagggatggacagttttacatgaggctgtgagcactggcgatcctgagatggtgtacacagtttccaacatcgagactaccacaacacatccatggcccttgaggagttcctgagctgctccaaaaattctcgaggctccggatt
tctatgtgcagatgaaatgggaattcaccagctgggtgcccttggtttctagaatatcccgaatgatgtctgtcgcacatcgaaaaagtgtgtgccaaactgcgcgtcgataacacattgctgggatttgaataacatgagctggataagaggaggcgtagttttatattaaaggagaagacaac
tggcgagtttaaggagtaacatgatgacaaagtgtgcaccaccgaacgttcgacctttccaagaatggagcgccctcactctggactgatgaagccaaaagcagggaagttagcggcggtcacaaagccctgtcattaacaccagctcgatactaaaaattgcttttga
aagggttttcagagttctcaagctaactctttgcaactgacagtgttagctga