

LOC402176-similar to 60S ribosomal protein L21 Gene

Reference: PTXBC057822

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms			
For research use only.			

GENE DESCRIPTION	similar to 60S ribosomal protein L21
GENE ID	402176
ORIGIN SPECIES	Human
ACCESSIONS	BC057822
VECTOR	pDONR223
NCBI SYMBOL	LOC402176
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggactgttccttgggtacaaatgcgaatctataagaaaggtgatatggtacacatcaacggaatggatactgttcaaaaagggaatgccccacaagtgttaccatggcaaaactgaaagtctacagtgttccccagcatgctgtggcattgtgtaaacaagggttaagattcttccaaga
gaattaatgtcatattgagcacatttaagcactgtaagagccgagatagcttctgaaatgtgtgaaggaaaatgacagaaaaagaaagccaaagagggaaggtacctgggttcaagtgaagcaccagcctgtcccccagagaaagcactttgtgaaagctgtgggaaggagc
ctgagctgctggaacctattccctgtgaattcatggcatcatag