

KRTAP26-1-keratin associated protein 26-1 Gene

Reference: PTXBC057825

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	keratin-associated protein 26-1; keratin associated protein 26-1		
For research use only.			

GENE DESCRIPTION	keratin associated protein 26-1
GENE ID	388818
ORIGIN SPECIES	Human
ACCESSIONS	BC057825
VECTOR	pDONR223
NCBI SYMBOL	KRTAP26-1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtcttgcctcaactactgtcgggaaactcaactcaggatctctcagaacctccgccatattctctacatccatcgacctctgccctaccagcgtgagctgtggagatgtcctctacttaccactagctctcaagaccatacctgggtcacagacaactgccaagagacctgcggtg
aaccaaccagctgccagccggtccactgtgagaccggcaatcttgaacctcttgcgggtttctcactgttactacgtgccagggccttgcgaaggaagcagttttctctctgtctttcttctcagttccttccagtagatctgttagaccacagagfatgtgtccagcggctgtctcca
ctgaggccgctgtcaatagttaccagcccataggagactgtgtgccaatgctatcgcccccaattctgttctccaagagttgccagccccaaacctctcactcttggatgccaacctcgagttgcttggcctatcgtcctcaaaagtcttccagctgtgtccagcagcctcagacctg
gggcctctgttcagtggttgccaacctctgacctatgtgttcagcagctgtgtccatcttctctgactgtga