



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

FAM3B-family with sequence similarity 3, member B Gene

Reference: PTXBC057829

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	protein FAM3B; C21orf11; C21orf76; ORF9; PRED44; D21M16SJHU19e; cytokine-like protein 2-21; pancreatic derived factor; family with sequence similarity 3 member B		

For research use only.

GENE DESCRIPTION	family with sequence similarity 3, member B
GENE ID	54097
ORIGIN SPECIES	Human
ACCESSIONS	BC057829
VECTOR	pDONR223
NCBI SYMBOL	FAM3B
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcccaattggctgggctctcaaggtgggtgtgtctctgcctctgtgtgcctggatctggggtaactctgtcgcagagatctatccagatgcacccctgtccagfctgcctatagatcctgcagatcggggagaggctgtgctctcaaaagtctcagctcccaaaaggcaa
aaatgtgacactggactcctggccactcagactcaggttactgctcaggggtgcagaagcaagtacgcacaaatctgctcttgaggataactactatgggagaaacagctcgggaaatgtgtccagaggataaataactgtcattgtcaactatgtactgggaatgtgcaga
cgaacacagatgtttgatgtatgaaggcgcataactctggacgcgatgacaagaagttatccagagctctcccaaaatccctgtctgtcgtgacactatgcagcgggaagcagactgaaataacgatgacagaagatgcatccaaagctgtcgaatgaagaatcgaacatgaaa
ttcagctctagctgtgtatgtatgcagcaaaaaggctgtgaactctggctgaaatctgcagcagaagaagaatgactgtctctgaagacaacagatattgtctggctgcagagatcagaagaagctgtcactatcccaaaagcgaagctgca

More info: Please contact us for any question: contact@proteogenix.fr