

DERL3-Der1-like domain family, member 3 Gene

Reference: PTXBC057830

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C22orf14; IZP6; LLN2; derlin-3; DERtrin 3; Der1-like domain family, member 3; degradation in endoplasmic reticulum protein 3; der1-like protein 3; derlin 3		

For research use only.

GENE DESCRIPTION	Der1-like domain family, member 3
GENE ID	91319
ORIGIN SPECIES	Human
ACCESSIONS	BC057830
VECTOR	pDONR223
NCBI SYMBOL	DERL3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcgtggcagggactagcggcaggttcctgacgtgcccgggtgacgcgggcttacaccgcagcctgtgtctcaccaccgcccgggtgcagctggagctcctcagcccttcaactctactcaaccgcacacctgtgtccggaagttccaggtctgaggctcgtcaccac
ttcctctcttcgggccccctgggaltcagcttcttcaacatgctctctgtgtccgctactgccgcatgctggaagagggtccttccgcccgcacggccgactctgtctcatgtttctcttcggggcgtccttatgacctgctgggactcctgggcagcctgttcttctggccaggc
cctcatggccatgctgtgtacgtgtgagccgcccagccctcgggtgagggtaacttcttccgctgctcacttccaggcaccgttctgcttggggcgtcatgggcttctcgtctgtctgggcaactccatcctcgtggacctgctgggaltgctgggcccatactactatc
ctggaggacgtcttcccaaccagcctggaggcaagagctcctgcagaccctggcttcttaagctgctcctggatgccctgcagaagaccctaattacctgccctcctgaggaaacagccaggaccatctgccaccccgagcagctga