

ZDHHC24-zinc finger, DHHC-type containing 24 Gene

Reference: PTXBC057833

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	DHHC-24; zinc finger DHHC domain-containing protein 24; zinc finger, DHHC domain containing 24; zinc finger DHHC-type containing 24		

For research use only.

GENE DESCRIPTION	zinc finger, DHHC-type containing 24
GENE ID	254359
ORIGIN SPECIES	Human
ACCESSIONS	BC057833
VECTOR	pDONR223
NCBI SYMBOL	ZDHHC24
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

```
atggggcagccctggggcctgggagcacggacgggcccgcgcagctgcctctctgctcaccgcgctgtgggcccggcctggcctggcttacgtgtgtgtctgtgtccggccgcgcctgggaccctggcccggccttcagctggcgtggcc
gccttcagctgctcaactctgctgggcaactgggctctctctgctcggatccagcatcctggcgtgatgtggccggccggctgtggccaggctggctfactgctaccaatgccaagccaggtgcgcccagcagcggacactgctctgctgcccgtctcatct
gcgtcgggaccaccactgccgcctgctggccgctgctgggcttcggcaactaccggccttctctgctcctgctctcatgccgcggcctgctctccactctctgtgctgctggccctgcaactgctggccctgctgcgagcccacacgccctccatggctgccctctctg
cttcctgctcatgtgtctcacaggcagagtgctctggcagatgcttgccttggccttgcgtgacggacacgtgcgtggcgggtgcgctgtgctgggctggcctgctcttccatgggatgctgctgctgcccagaccacatgggagtggtcgggcccagcactctatgacc
tgggtccctgccacaacctgcaggcagccctgggcccctgctgctgctgcccctctgctcctcccaattgctggggatgggatcacctccagaccacagcagatgggacacacagcctctga
```