

TTC32-tetratricopeptide repeat domain 32 Gene

Reference: PTXBC057850

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	tetratricopeptide repeat protein 32; TPR repeat protein 32; tetratricopeptide repeat domain 32		

For research use only.

GENE DESCRIPTION	tetratricopeptide repeat domain 32
GENE ID	130502
ORIGIN SPECIES	Human
ACCESSIONS	BC057850
VECTOR	pDONR223
NCBI SYMBOL	TTC32
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggaaggacagcggcaagaagccacgcaaccctaactcgcggcgcctcattcaacaatggagagtagcgggaggcggaggcactgtactcgccttacatcgcgggtgcgttgcgcggcctccagcgacgagagtcgggagcaaatgcagccctgaggatttgctact
gcataatacaacagggggcaaatcaagtactcagggttgattttatgaagccatggatgactacacatctgcatagaagtccaaccaattttgaagtccatattacaacagagggttgatactgtataggctgggatatgttgatgatccttgggaagattcaagaaggctcttagactaaat
cctggattcaagatgctactttgagcttaaacagactattctagacaagaagaaaaacaaagaagaatgttgcaaaaaattattga