

LYSMD3-LysM, putative peptidoglycan-binding, domain containing 3 Gene

Reference: PTXBC058027

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	lysM and putative peptidoglycan-binding domain-containing protein 3; LysM, putative peptidoglycan-binding, domain containing 3; LysM domain containing 3		

For research use only.

GENE DESCRIPTION	LysM, putative peptidoglycan-binding, domain containing 3
GENE ID	116068
ORIGIN SPECIES	Human
ACCESSIONS	BC058027
VECTOR	pDONR223
NCBI SYMBOL	LYSMD3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcaggaggagcagcaatcgtagtttctctccaggagttcagtcagtggtcaagtacatgcatttgaaattgttcagacagtgatatttgaggaggatgctgaagtgtatgaactcgcagcagagagaaagagaaagtcggaagaagatcatcaagagatagacttgacgaca
ttatagtattacaagaagatatacaagaaggagatacaattaaatgcaatagcccttcagtcactgtgtacgggtctatcaaaattccagtaaaaaagttcagttccttgaccgaaacactttgtctccaaaaggagacagactcagctcattcctgttcaatactctccgaacaacaggaaattt
tgccagctaa