



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

HP-haptoglobin Gene

Reference: PTXBC058031

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	HP2ALPHA2; HPA1S; binding peptide; haptoglobin alpha(1S)-beta; haptoglobin alpha(2FS)-beta; haptoglobin, alpha polypeptide; haptoglobin, beta polypeptide; zonulin		

For research use only.

GENE DESCRIPTION	haptoglobin
GENE ID	3240
ORIGIN SPECIES	Human
ACCESSIONS	BC058031
VECTOR	pDONR223
NCBI SYMBOL	HP
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggtttccccaataaactaccacaggtgccacgcgtatcaatgacaatgctgctgcaccagcgctaaaaatctctctgcgaacattcagaaaaatgcaacagcgaagacatgcccctactttfaacctctatgtgtgggaaaaagacgctgtagagatgagaagggtttctacacc
taactactccaggtagatattgggctcatcaaacccaacagaagggtgtctttaaagagagagtgatgccatctgcctactctcaaggattatgcagaagttagggcggtgtgtgggtgtatttctgcgtgtggggcgaaatgcacaaatgaaatctgaccatctgaagtatgctatgctgct
gtggctgcacacacaaatgcataaggcattatgaagacagacagatccccgaagaagaagacacagaagacgcctgtagggggtgcagccatcatgaatgaacacacatctgtgcgtatctgaagtacagaagaacacacctgctatgcgcatgcggcagtcgctttgcccttcac
gacctgcagcgaacacatctctaccacacatctaaactctaaactctgaaagaaatctgtctgcctgcgaatgtatgaaggtgcctctcatccacgaactgcgttcagaaagacacatactgctgaataactaa