



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

C14orf118-chromosome 14 open reading frame 118 Gene

Reference: PTXBC058032

Product type:

DNA

Storage Condition:

-20°C

Size/Quantity:

2 ug

Delivery Condition:

+4°C, Blue Ice

Synonyms

C14orf118; G patch domain-containing protein 2-like; G-patch domain containing 2 like

For research use only.

GENE DESCRIPTION	chromosome 14 open reading frame 118
GENE ID	55668
ORIGIN SPECIES	Human
ACCESSIONS	BC058032
VECTOR	pDONR223
NCBI SYMBOL	C14orf118
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr