

RASL10A-RAS-like, family 10, member A Gene

Reference: PTXBC058077

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	ras-like protein family member 10A; ras-like protein RRP22; ras-related protein on chromosome 22; RAS like family 10 member A		

For research use only.

GENE DESCRIPTION	RAS-like, family 10, member A
GENE ID	10633
ORIGIN SPECIES	Human
ACCESSIONS	BC058077
VECTOR	pDONR223
NCBI SYMBOL	RASL10A
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgggggtagcctcgggtggcgtttaggcgccccgggctgggcaagacggccatcatccagttctgttcggtgactacccgagcgccacggccacggagcggccgctctaccgaccgctgtctcgcagggcgccgtctacgactgagcatccgga
cggcgacgtcgtgccccggctcgagccccgggggtcggagagtgccagacgctaaggactggagcttgaggacacggagccttcgtctgtctacgacatctgcagcccgacagtttcgactacgtgaaggccctgcggcagcgcatcgaggaccgagccggc
ggcgcgccccgaagcgccatcctcgtgtgtaggcaacaagcgggacaggcagcggctgcgttcggaccgcgccgctgtgccccttagtgccgagggctggcgctgcggctacctcagtgctccgccaagtacaactggcagctgtcgtctctccgagctgtg
cgtgcgtctgtgctgcgcgcgcctgcacaccggccctgcctgcaggggcgctgcatcccgcgctgcagcctcatgtga