

AGR3-anterior gradient homolog 3 (*Xenopus laevis*) Gene

Reference: PTXBC058284

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	AG-3; AG3; BCMP11; HAG3; PDIA18; hAG-3; anterior gradient protein 3; anterior gradient homolog 3; anterior gradient protein 3 homolog; breast cancer membrane protein 11; protein disulfide isomerase family A, member 18; anterior gradient 3, protein disulph		

For research use only.

GENE DESCRIPTION	anterior gradient homolog 3 (<i>Xenopus laevis</i>)
GENE ID	155465
ORIGIN SPECIES	Human
ACCESSIONS	BC058284
VECTOR	pDONR223
NCBI SYMBOL	AGR3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgatgtacactcagcttgggtctctgctcttactctgacagtttctccaacctgccattgcaataaaaaaggaaaaggccctctcagacactctcaagaggatggggagatgacatcacttgggtacaaacttatgaagaaggtctctttatgtcaaaaaagtaagaagccattaa
tgggtattcatcacctggaggattgtcaatactctcaagcactaaagaaagtatttgcacaaatgaagaatacaagaataatggctcagaataagttcatcatgctaaaccttatgcatgaaccactgataagaattatcactgatgggcaatattgctcctagaatcatgtttgtagaccctcttt
aacagtttagagctgacatagctggaagatacttaacagattgtacacatatgagcctcgggatttaccctattgatagaacaatgaagaagcattaagacttattcagtcagagctataa