

LOC204010-ribosomal protein SA pseudogene Gene

Reference: PTXBC058822

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms			
For research use only.			

GENE DESCRIPTION	ribosomal protein SA pseudogene
GENE ID	204010
ORIGIN SPECIES	Human
ACCESSIONS	BC058822
VECTOR	pDONR223
NCBI SYMBOL	LOC204010
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgctgaagtttgcctgacctggagccactccaattgctggccacttcaactctggaaccttctaaccagatccaggcagcctactgggagccacggcttctgtgtttctgacccagggtgaccaccagtcctcacagaggcatctgtgttaaccacctgccattgctttatgca
acacagattctcctctgtgccattggacattgccatcacatgcaacaacaaggagctccctcagtggtgagatgtgttagatgctggcctgggaagtctgcacatgctggcaccatttctgcaacacccatggaggtaatgtctaatctctactctacagagatcctgaagagactg
aaaaagaagagcagggtgctgtgaaaaaggctgtgaccaagaaggaaatttcagggttaattggactgttccagctcctgagttcactgttactcagcctgaagttgcagactggctgaaaggcatgagggtgcccctctgtgcctattcagcagtttctactgaagactggaaggctcagcctgc
catggaagactggtccacagctccactgctcagccactga