

GOLT1A-golgi transport 1 homolog A (*S. cerevisiae*) Gene

Reference: PTXBC058832

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CGI-141; GOT1; HMFN1187; YMR292W; hGOT1b; vesicle transport protein GOT1A; golgi transport 1A		

For research use only.

GENE DESCRIPTION	golgi transport 1 homolog A (<i>S. cerevisiae</i>)
GENE ID	127845
ORIGIN SPECIES	Human
ACCESSIONS	BC058832
VECTOR	pDONR223
NCBI SYMBOL	GOLT1A
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgatctccatcaccgaatggcagaagattggtgtgggacaccggtttcggcatcttctcatctctttggaacactcctgtactttgattcgtgctcctggcctttggaaacctgctgttctgacgggctgtccctcatcattggcctgaggagacctttgttctttccaacggcaca
aactcaagggaaccagcttctcctgggggtgtgtatcgtgctcctacgtggccctcctcgcatgttcttggaacctacggattcttcagctctttaaggcctttttccctgtcgccttcggcttctgggcaatgtctgcaacatccccctctgggtgcgctgttccggagactcaa
ggcactagctcgatggtctga