

WDR3-WD repeat domain 3 Gene

Reference: PTXBC058836

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	DIP2; UTP12; WD repeat-containing protein 3; dJ776P7.2 (WD repeat domain 3); WD repeat domain 3		

For research use only.

GENE DESCRIPTION	WD repeat domain 3
GENE ID	10885
ORIGIN SPECIES	Human
ACCESSIONS	BC058836
VECTOR	pDONR223
NCBI SYMBOL	WDR3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggggctcaccaagcagctacgtatgttgctagtgcggtcttggcggtatcgagccaaaaggtctagtagccaagcatcaagttccttgtgttcttcgccactgaattctga