

ADAMTS12-ADAM metalloproteinase with thrombospondin type 1 motif, 12 Gene

Reference: PTXBC058841

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	PRO4389; A disintegrin and metalloproteinase with thrombospondin motifs 12; ADAM-TS 12; ADAM-TS12; ADAMTS-12; a disintegrin-like and metalloprotease (reprolysin type) with thrombospondin type 1 motif, 12; ADAM metalloproteinase with thrombospondin type 1 m		

For research use only.

GENE DESCRIPTION	ADAM metalloproteinase with thrombospondin type 1 motif, 12
GENE ID	81792
ORIGIN SPECIES	Human
ACCESSIONS	BC058841
VECTOR	pDONR223
NCBI SYMBOL	ADAMTS12
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgccatgtgccagaggagctggctgcaaaccttccgtgggtcagctccttaacttggggcgcttgcctatgggagacagcctcagcaggccgggttcgctcccgacaggaggaagagcatttatcaaggcctgccagaataaccagtggtgggtccagtcagtag
atgccagtgggcatttttgcataatggcttgcactatcccatcacgagcagcagggaagagagatttggatggctcagaggactgggtgactacagaatttctcagaggagaaggacgtgttttaacttgacgggtcaatcaaggatttcttccatagctacatcatggagaagagat
atgggaacctctccatgttaagatgatggcttctcagccccctctgccatctcagtcgacaggttctcagcagggcaccagagttgggacggcagccctcagtgccctgcatggactgactggattttccaactaccatggagacttttcaattgaacccgtgaagaagcatccact
ggttgaggagggtaccaccgcacatcgtttacaggaggcagaaagtccagaaccaaggagccaactgtgattaaagggtattgtgactcacatgtcctcctgggttgaagaatctgtttgtctttgttag