



Reference: PTXBC058844

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	HLUGP4; LUGP4; UDPGT; UDPGT 1-9; UGT-1I; UGT1-09; UGT1-9; UGT1.9; UGT1A9S; UGT1AI; UGT1I; UDP-glucuronosyltransferase 1-9; UDP glucuronosyltransferase 1 family, polypeptide A9; UDP glycosyltransferase 1 family, polypeptide A9; UDP-glucuronosyltransferase		

For research use only.

GENE DESCRIPTION	UDP glucuronosyltransferase 1 family, polypeptide A9
GENE ID	54600
ORIGIN SPECIES	Human
ACCESSIONS	BC058844
VECTOR	pDONR223
NCBI SYMBOL	UGT1A9
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]