

**SULT1C4-sulfotransferase family, cytosolic, 1C, member 4 Gene**

Reference: PTXBC058861

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	SULT1C; SULT1C2; sulfotransferase 1C4; ST1C4; SULT1C#2; sulfotransferase 1C2; sulfotransferase family, cytosolic, 1C, member 2; sulfotransferase family, cytosolic, 1C, member 4; sulfotransferase family, cytosolic, 1C, member C2; sulfotransferase family 1C		

**For research use only.**

---

<b>GENE DESCRIPTION</b>	sulfotransferase family, cytosolic, 1C, member 4
<b>GENE ID</b>	27233
<b>ORIGIN SPECIES</b>	Human
<b>ACCESSIONS</b>	BC058861
<b>VECTOR</b>	pDONR223
<b>NCBI SYMBOL</b>	SULT1C4
<b>LENGTH</b>	ProteoGenix bp
<b>SEQUENCE PRIMERS</b>	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

**ORF (Only for reference, product provided without STOP codon )**

atggccttacacgacatggaggatttacattgatggaacaaagcgttaagtgtcaactacgtgaagggaattctcaaccgacagacacctgtgacatctggataagatctggaactccaagccaagcctgatgacctgcttattctacctatcctaaagcaggaacaacatggactca  
ggagatagtggaattaatacaaaatgaaggatgagagaaaagtaaacggcaccgactcatcaacgatttccctcctcgaaatgaaaatcccatccttaggatctgtgagtataaatag