



Reference: PTXBC058900

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	FRA; rab GTPase-binding effector protein 2; Rabaptin-5 beta; rabaptin-5beta; rabaptin, RAB GTPase binding effector protein 2		

For research use only.

GENE DESCRIPTION	rabaptin, RAB GTPase binding effector protein 2
GENE ID	79874
ORIGIN SPECIES	Human
ACCESSIONS	BC058900
VECTOR	pDONR223
NCBI SYMBOL	RABEP2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]