

## SFRP4-secreted frizzled-related protein 4 Gene

Reference: PTXBC058911

|                |                                                                                                                                                                                                                    |                     |                |
|----------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------|----------------|
| Product type:  | DNA                                                                                                                                                                                                                | Storage Condition:  | -20°C          |
| Size/Quantity: | 2 ug                                                                                                                                                                                                               | Delivery Condition: | +4°C, Blue Ice |
| Synonyms       | FRP-4; FRPHE; PYL; sFRP-4; secreted frizzled-related protein 4; frizzled protein, human endometrium; secreted frizzled-related protein 4; secreted frizzled-related protein 4; secreted frizzled related protein 4 |                     |                |

**For research use only.**

---

|                         |                                                                                                   |
|-------------------------|---------------------------------------------------------------------------------------------------|
| <b>GENE DESCRIPTION</b> | secreted frizzled-related protein 4                                                               |
| <b>GENE ID</b>          | 6424                                                                                              |
| <b>ORIGIN SPECIES</b>   | Human                                                                                             |
| <b>ACCESSIONS</b>       | BC058911                                                                                          |
| <b>VECTOR</b>           | pDONR223                                                                                          |
| <b>NCBI SYMBOL</b>      | SFRP4                                                                                             |
| <b>LENGTH</b>           | ProteoGenix bp                                                                                    |
| <b>SEQUENCE PRIMERS</b> | Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT<br>Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG |

### ORF (Only for reference, product provided without STOP codon )

atgttctctccatctagtggcgtgtgctgtgctgcaactggcgctggcggtgcgcggcgccctgcagggcggtgcgcacatccctatgtgccggcacatgccctggaacatcacgcggatgcccaaccacctgcaccacagcacgaggagaacgccatcctggccatcagc  
agtagcaggagctggtggacgtgaactgcagcgccgtgctgcgtctctctctctgtgcatgtacgcgccattgcacccctggagttcctgcagcaccctatcaagccgtgcaagtcgggtgtgccaaagcgcgccgcgacgactgcgagccctcatgaagatgtacaaccacagctggc  
ccgaaagcctggcctgcgacgagctgctgtatgacgtggcggtgtgcatctgcctgaagccatcgtcacggacctcccgaggatgttaagtggatgacatcacaccagacatgatgtacaggaagccctctgtgtgactgtaaacgcctaaagcccgatcgggtgcaagt  
taaaaaggtgaagccaactttggcaactatctcagcaaaaactacagctatgttattcatgccaaaataaaagctgtgcagaggagtggtgcgaatgaggtcacacacggtgtgtggtgttaaaagagatcttcaagtcctcatcccatccctcgaactcaagtcctcattacaattct  
cttgcagtgccacacatcctgccccatcaagatgttctcatcatgtgttacgagtggtgtcaaggtatgctcttgaattgcttagtgaataatggagagatcagcttagtaaaagatccatacagtggaagagaggtgtcaggaacagcgagaaacaggttcaggacaagaaga  
aaacagccggcgccacagctgtagtaatcccccaaacaaagggaagcctctgctcccaaccagccagtcgaagaagaacattaaactaggagtgccgaagagaaacaaacccgaaaagagtggtga