

UBL4B-ubiquitin-like 4B Gene

Reference: PTXBC058929

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	ubiquitin-like protein 4B; ubiquitin like 4B		
For research use only.			

GENE DESCRIPTION	ubiquitin-like 4B
GENE ID	164153
ORIGIN SPECIES	Human
ACCESSIONS	BC058929
VECTOR	pDONR223
NCBI SYMBOL	UBL4B
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgttcctcacagtcaagctgctcctggccagagatgcagctgaaggtgtcaggccaagagatgtagccacgtgaagagactgggtccaggcggctgaaggtgcctgaggagcagcagcagctctttccgtggccagctcctggaggatgacaagcacctcttgactactgc
attgggccaatgcctctatcaatgtcatcatgcagccctggagaagatggcgctaaaggaggccaccagccgagaccagccctgtggcaccagctgggactggctctagctaaacacttgaaccacagatgccaaggccgtgctgagctgctaaggcaggagcagagg
agcgctgcagaagataagcctggagcacctggagcagctggccagctacctcctggcagaggagcctacgtggagccagctggagagaggagcttgaggcgaaggcacggcctcagagctcctgtgacatggaggagaaggaggagcagcagctgatcagtaa