



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

IKIP-IKK interacting protein Gene

Reference: PTXBC058933

Product type:

DNA

Storage Condition:

-20°C

Size/Quantity:

2 ug

Delivery Condition:

+4°C, Blue Ice

Synonyms

IKIP; inhibitor of nuclear factor kappa-B kinase-interacting protein; I kappa B kinase interacting protein; IKK interacting protein; i kappa-B kinase-interacting protein; IKBKB interacting protein

For research use only.

GENE DESCRIPTION	IKK interacting protein
GENE ID	121457
ORIGIN SPECIES	Human
ACCESSIONS	BC058933
VECTOR	pDONR223
NCBI SYMBOL	IKIP
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr

Page 1/1