



Reference: PTXBC059376

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	gamma-aminobutyric acid receptor subunit epsilon; GABA(A) receptor subunit epsilon; GABA(A) receptor, epsilon; gamma-aminobutyric acid (GABA) A receptor, epsilon; gamma-aminobutyric acid type A receptor epsilon subunit		

For research use only.

GENE DESCRIPTION	gamma-aminobutyric acid (GABA) A receptor, epsilon
GENE ID	2564
ORIGIN SPECIES	Human
ACCESSIONS	BC059376
VECTOR	pDONR223
NCBI SYMBOL	GABRE
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr