



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

FAM108C1-family with sequence similarity 108, member C1 Gene

Reference: PTXBC059401

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	abhydrolase domain-containing protein FAM108C1; FAM108C1; protein ABHD17C; abhydrolase domain-containing protein 17C; alpha/beta hydrolase domain-containing protein 17C; family with sequence similarity 108, member C1; abhydrolase domain containing 17C		

For research use only.

GENE DESCRIPTION	family with sequence similarity 108, member C1
GENE ID	58489
ORIGIN SPECIES	Human
ACCESSIONS	BC059401
VECTOR	pDONR223
NCBI SYMBOL	FAM108C1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr