



Reference: PTXBC059402

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	round spermatid basic protein 1-like protein; round spermatid basic protein 1 like		
For research use only.			

GENE DESCRIPTION	round spermatid basic protein 1-like
GENE ID	222194
ORIGIN SPECIES	Human
ACCESSIONS	BC059402
VECTOR	pDONR223
NCBI SYMBOL	RSBN1L
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr