

OVOL1-ovo-like 1(Drosophila) Gene

Reference: PTXBC059408

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	HOVO1; ovo homolog-like 1; ovo like zinc finger 1; ovo-like 1(Drosophila); ovo like transcriptional repressor 1		

For research use only.

GENE DESCRIPTION	ovo-like 1(Drosophila)
GENE ID	5017
ORIGIN SPECIES	Human
ACCESSIONS	BC059408
VECTOR	pDONR223
NCBI SYMBOL	OVOL1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgccccgcgcttctggtgaagaagccgtgcgtctccacgtgcaagaggaaactggagcagctccccgacgaggagcggcgagatctacgtgccagtcagcctgggtcttgcaccacagccctaccgggagccgggaacctctgtggcgaacccctctgcccgcgtg
gcttgaacatgagccttcgagactctagctacagcatggccccgggcccctgtgtgtggcccagctgccccttgaagacatgggccacttgacagacccccagagcagagaccatggcttctgcgcaccaagatgaaggfagacctggggacagfcccagtgagacctgttcac
ctgccgtgtctgccagaaggcccttcacctaccagcgcgtgctgaaccgcccacatgaagtgtcacaaacgacgtcaagaggcacctctgcacgtactgcgggaagggtctcaatgacaccttcgacctcaagagacacgtccgaactcacactggcgtgcggccctacaagtgcagcctgt
gtgacaaggccttcacgcagcgtgctctctggagctcacctcaagaagatccatggtgtgcagcagaagtacgcgtacaaggagcggcgcccaagctgtactgtgtgaggagtgccgctgcacatctgagagccaggaggccacgtctgcacctgaaggagcaccacctg
acagcccgtctgctgcgaagacctccaagaaggtggccgtggcactacagaacactgtcacttcctgctgcagggcagccccacctgtga