

GKN1-gastrokine 1 Gene

Reference: PTXBC059778

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	AMP18; BRICD1; CA11; FOV; foveolin; gastrokine-1; 18 kDa antrum mucosa protein; AMP-18; BRICHOS domain containing 1; gastrokine 1		

For research use only.

GENE DESCRIPTION	gastrokine 1
GENE ID	56287
ORIGIN SPECIES	Human
ACCESSIONS	BC059778
VECTOR	pDONR223
NCBI SYMBOL	GKN1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaagttcacaattgtcttctgctggactcttggagctcttctagctcctgcctagctactataatcaacgtcaatgatgacaacaacaatgtggaagtgggcagcagtcagtgagtgcaacaatgaacacaatgtggccaatgttgacaataacaacggatgggactcctggaattcc
atctgggattatggaaatggctttctgcaaccagactcttcaaaaagaagacatgctgtgcacaaatgaacaaggagtcacccctccatcaatcccttgatgcactgggtcaaggaaaagaagctcagggttaagggaaccaggaggaccacctccaaaggccctgatgtactcagt
caaccacaacaagtcgatgacctgagcaagtcgggaaaaacattgcaaacatgtgtcgtgggattccaacatacatggctgaggagatgcaaggagcaagccctgttttttactcagggaacgtgctacacgaccagtgactatgattgtggacattcctctgtggagacacgggtgga
gaactaa