

LGALS2-lectin, galactoside-binding, soluble, 2 Gene

Reference: PTXBC059782

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	HL14; galectin-2; S-Lac lectin 2; beta-galactoside-binding lectin L-14-II; gal-2; lactose-binding lectin 2; lectin, galactoside-binding, soluble, 2; galectin 2		

For research use only.

GENE DESCRIPTION	lectin, galactoside-binding, soluble, 2
GENE ID	3957
ORIGIN SPECIES	Human
ACCESSIONS	BC059782
VECTOR	pDONR223
NCBI SYMBOL	LGALS2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgacgggggaacttgaggtaagaacatggacatgaagccgggtcaaccctgaagatcacaggcagcatgccgatggcactgatggctttgtaattaatctgggaggggacagacaagctgaacctgcattcaaccctcgttcagcgaatccaccattgtctgcaactcattgg
acggcagcaactggggcaagaacaacgggaagatcacctgtgttcagcccagggtcagagggtcaagttcacagtgacctttgagagtgcacaaattcaagggtgaagctgccagatgggcagcagctgacttttccaacaggctgggtcacagccacctgagctacctgagcgttaag
gggcgggttcaacatgtcctcttcaagttaaaagaataa