

CSAG1-chondrosarcoma associated gene 1 Gene

Reference: PTXBC059947

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CT24.1; cancer/testis antigen 24.1; cancer/testis antigen CSAGE; cancer/testis antigen family 24, member 1; chondrosarcoma associated gene 1		

For research use only.

GENE DESCRIPTION	chondrosarcoma associated gene 1
GENE ID	158511
ORIGIN SPECIES	Human
ACCESSIONS	BC059947
VECTOR	pDONR223
NCBI SYMBOL	CSAG1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtcggcgactacagcctgtcgtgcctcactgtcctgggggaagctcggggagaccaggtggactggagtagactgtacagagacactggtctggtgaagatgtccaggaaaccacgagcctcagccattttccaacaaccacccatcaacaccaaagaggtccaagaca
accagaagggaagggacccgtcaagggaagttccaggaaacaaaggtctccctaa