

RPL12-ribosomal protein L12 Gene

Reference: PTXBC059950

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	L12; 60S ribosomal protein L12; ribosomal protein L12		

For research use only.

GENE DESCRIPTION	ribosomal protein L12
GENE ID	6136
ORIGIN SPECIES	Human
ACCESSIONS	BC059950
VECTOR	pDONR223
NCBI SYMBOL	RPL12
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgccgcgaagtgcacccaacgagatcaaatcggtatgaggtgacccggaggtgaagtcggtgccacttctgccctggccccaagatcgcccccgggtctgattgaggtgggtgcctctgcctctgccctgatcatcaagccctcaaggaaaccacagagacagaaag
aaacagaaaaacattaaacacagtggaatatactttgatgagattgtcaacattgctgcacagatgcggcacgcgacgttagccagagaactcttgaaccattaaagagatcctggggactgccagtcagtggtgtaattgtgatggccgccatcctcatgacatcatcgtgaca
tcaacagtggtgctgtggaatgccagccagttaa