

C19orf33-chromosome 19 open reading frame 33 Gene

Reference: PTXBC060319

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	H2RSP; IMUP-1; IMUP-2; immortalization up-regulated protein; HAI-2 related small protein; hepatocyte growth factor activator inhibitor type 2-related small protein; immortalization-upregulated protein; chromosome 19 open reading frame 33		

For research use only.

GENE DESCRIPTION	chromosome 19 open reading frame 33
GENE ID	64073
ORIGIN SPECIES	Human
ACCESSIONS	BC060319
VECTOR	pDONR223
NCBI SYMBOL	C19orf33
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggagttcgacctgggagcagccctggagcccacctcccagaagcccggtgtggggcgggccacgggggagatcccaagctcagtcgcccaaaagttcaggccggtcggagcaggggcaggtccgggtccaaaggcctccaactcagggggctgggtaagggcgccg
cctcactgcgcacctccatccagcaaggacaccacagctctccgactccagcagcagctccagcgattcggacacggatgtga